

Product Datasheet

Anti-h CEA 5909 SP-5 100428

Product Name	Anti-h CEA 5909 SP-5
Catalog Number	100428
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Note	In Bhayana et al. (1989) authors made an epitope mapping for 11 monoclonal antibodies from Medix Biochemica and they found that the antibodies bound to 5 different epitope groups (A to F) and 5909 was found to bind to epitope group C.
Isotype	IgG1
Clonality	Monoclonal
Epitope	Epitope group C as described in Bhayana et al. (1989)
Purity	≥ 95 %
Affinity constant	$KA = 3 \times 10^{10}$ 1/M
Associated Products	Native CEA antigen, Lee Biosolutions 151-09 and 151-10
Buffer	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.0-7.1
Cross Reactivity	Does not recognize human NCA (non-specific cross-reacting antigen), NCA-2 or BGP-1 (biliary glycoprotein-1)
Specificity	Antibody recognizes human carcinoembryonic antigen
Shelf Life	36 months
References	Bhayana, V. and Diamandis, E.P. (1989) A double monoclonal time-resolved immunofluorometric assay of carcinoembryonic antigen in serum. Clin. Biochem., 22:433-438