

Product Datasheet

Anti-Human IL-29 Monoclonal Antibody, Azide Free Clone B-Z29 855.880

Product Name	Anti-Human IL-29 Monoclonal Antibody, Azide Free Clone B-Z29
Catalog Number	855.880
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant human IL-29
Alternative Names	Interferon lambda 1 ZCYTO21
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes both natural and recombinant human IL-29
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells