

Product Datasheet

Anti-Human CD196 Monoclonal Antibody, Azide Free Clone B-R35 855.910

Product Name	Anti-Human CD196 Monoclonal Antibody, Azide Free Clone B-R35
Catalog Number	855.910
Description	Mouse monoclonal antibody. Source: SP2/O x Balb/c B cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	huCCR6 NIH-3T3 transfected cells
Alternative Names	CCR6
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	CD196
Hybridoma	SP2/O x Balb/c B cells