

Product Datasheet

In vivo IgG2b Isotype Control, Azide Free Clone B-F33 855.940.005

Product Name	In vivo IgG2b Isotype Control, Azide Free Clone B-F33
Catalog Number	855.940.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Control
Immunogen	Chemical molecule
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG2b Lambda light chain
Clonality	Monoclonal
Specificity	Isotype control for in-vivo experimentation
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells