

Product Datasheet

In vivo IgG1 Isotype Control, Azide Free Clone B-D38 855.950.005

Product Name	In vivo IgG1 Isotype Control, Azide Free Clone B-D38
Catalog Number	855.950.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Control
Immunogen	Chemical molecule
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Isotype control for in-vivo experimentation
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells