
Product Datasheet

Human IL-6 ELISpot Set 856.021

Product Name	Human IL-6 ELISpot Set
Catalog Number	856.021
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent. Spots read visually or with reader.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes natural human IL-6
References	Kalogerakou, F. et al.,Hippokratia,2008; 12(4): 230-5. - Pubmed link