

Product Datasheet

Human IL-10 ELISpot Set 856.031.015

Product Name	Human IL-10 ELISpot Set
Catalog Number	856.031.015
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, and blocking reagent.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes natural human IL-10