

## Product Datasheet

### Human TNF alpha ELISpot Set 856.041

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Human TNF alpha ELISpot Set  |
| <b>Catalog Number</b>      | 856.041  |
| <b>Description</b>         | Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent. |
| <b>Tested Applications</b> | ELISpot  |
| <b>Alternative Names</b>   | TNFa   |
| <b>Brand</b>               | Diaclone   |
| <b>Storage</b>             | Store at +4°C  |
| <b>Note</b>                | Spots read visually or with reader   |
| <b>Cross Reactivity</b>    | No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with TNFR I and TNFR II.                      |
| <b>Specificity</b>         | Recognizes natural human TNF-a   |