

## Product Datasheet

### Human TNF alpha ELISpot Set 856.041.015

<b>Product Name</b>	Human TNF alpha ELISpot Set
<b>Catalog Number</b>	856.041.015
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, and blocking reagent.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	TNFa
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with TNFR I and TNFR II.
<b>Specificity</b>	Recognizes natural human TNF-a