

Product Datasheet

Human TNF alpha ELISpot Set 856.041.020

Product Name	Human TNF alpha ELISpot Set
Catalog Number	856.041.020
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, and blocking reagent.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with TNFR I and TNFR II.
Specificity	Recognizes natural human TNF-a