

Product Datasheet

Human IL-17A ELISpot Set 856.151

Product Name	Human IL-17A ELISpot Set
Catalog Number	856.151
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent. Spots read visually or with reader.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IL-17A
References	Rinaldi, M. et al., Thorax, 2012; 10.1136/thoraxjnl-2011-200690 - Pubmed link