

Product Datasheet

Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12 857.000

Product Name	Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12
Catalog Number	857.000
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC, IP
Immunogen	KG-1 cell line
Alternative Names	PRR2 HVEB
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises Nectin-2, the polio retrovirus receptor-2 (PRR2)
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells