

Product Datasheet

Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12 857.000.005

Product Name	Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12
Catalog Number	857.000.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IP, FC
Immunogen	KG-1 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Nectin-2, the polio retrovirus receptor-2 (PRR2)
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells