

Product Datasheet

Anti-Human Neurotensin R1 Monoclonal Antibody, Azide Free Clone B-N6 857.040

Product Name	Anti-Human Neurotensin R1 Monoclonal Antibody, Azide Free Clone B-N6
Catalog Number	857.040
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC, Functional assay
Immunogen	NTR transfected CHO cell line
Alternative Names	Neurotensin R1
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. DO NOT FREEZE. DO NOT CONCENTRATE (risk of precipitation).
Isotype	Mouse IgM Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises the membrane protein
Biological Activity	Inhibits neurotensin binding to its receptor
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells