

Product Datasheet

Anti-Biotin Monoclonal Antibody, Azide Free Clone B-E24 857.060.005

Product Name	Anti-Biotin Monoclonal Antibody, Azide Free Clone B-E24
Catalog Number	857.060.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, WB, IHC
Immunogen	KLH conjugated Biotin
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises biotin
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells