

Product Datasheet

IgG1 Isotype Control, Azide Free Clone B-Z1 857.070.005

Product Name	IgG1 Isotype Control, Azide Free Clone B-Z1
Catalog Number	857.070.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC, FC
Immunogen	Rat monocytes and macrophages
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Rat monocytes and macrophages
Specificity	Negative control for biological test
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells