

Product Datasheet

Anti-Human IL1RAP Monoclonal Antibody, Azide Free Clone B-L43 857.090

Product Name	Anti-Human IL1RAP Monoclonal Antibody, Azide Free Clone B-L43
Catalog Number	857.090
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant human IL1RAP
Alternative Names	IL-1 RAcP (Interleukin-1 Receptor Accessory Protein) IL1 R3 C3orf13
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes natural and recombinant human IL1RAP
References	Trad, R. et al., J Immunother Cancer. 2022; 10(7): e004222.
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells