

Product Datasheet

Anti-Human IL1RAP Monoclonal Antibody, Unconjugated Clone B-R58 857.113.020

Product Name	Anti-Human IL1RAP Monoclonal Antibody, Unconjugated Clone B-R58
Catalog Number	857.113.020
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human IL1RAP
Alternative Names	IL-1 RAcP (Interleukin-1 Receptor Accessory Protein) IL1 R3 C3orf13
Brand	Diaclone
Form/Appearance	Liquid
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognizes natural and recombinant human IL1RAP
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells