

Product Datasheet

Anti-SARS-CoV-2 Spike S1 Protein Azide Free Clone CR3022 857.140.005

Product Name	Anti-SARS-CoV-2 Spike S1 Protein Azide Free Clone CR3022
Catalog Number	857.140.005
Description	Mouse monoclonal
Species/Host	Human
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV.
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Human IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from a phage display library
Specificity	S1 domain of the SARS-CoV Spike protein (aa318-510) as well as SARS-CoV-2 (COVID-19) Spike protein