

Product Datasheet

Anti-SARS-CoV-2 Nucleoprotein Azide Free Clone CR3018 857.150

Product Name	Anti-SARS-CoV-2 Nucleoprotein Azide Free Clone CR3018
Catalog Number	857.150
Description	Source: Recombinant antibody from a phage display library
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV
Alternative Names	Nucleoprotein SARS-CoV nucleoprotein Nucleocapsid protein
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Human IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from a phage display library
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Nucleoprotein of SARS-CoV-2 (aa 11-19) which corresponds to the sequence RSAPRITFG of the N protein of SARS-CoV