

Product Datasheet

Anti-SARS-CoV-2 Spike S1 Protein Azide Free Clone CR3022 (IgA) 857.160.005

Product Name Anti-SARS-CoV-2 Spike S1 Protein Azide Free Clone CR3022 (IgA)

Catalog Number 857.160.005

Description Mouse monoclonal

Species/Host Human

Conjugation Azide free

Tested Applications ELISA

Immunogen Generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV.

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months.

Isotype Human IgA Kappa light chain

Clonality Monoclonal

Source Recombinant antibody from a phage display library

Specificity S1 domain of the SARS-CoV Spike protein (aa318-510) as well as SARS-CoV-2 (COVID-19) Spike protein