

Product Datasheet

Anti-Human FGF2 Monoclonal Antibody,Azide Free Clone B-T42 857.190

Product Name	Anti-Human FGF2 Monoclonal Antibody,Azide Free Clone B-T42
Catalog Number	857.190
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human FGF basic
Alternative Names	HBGF-2 FGF basic Fibroblast Growth Factor Basic
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes natural and recombinant human FGF2
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells