

## Product Datasheet

### Anti-Human HLA Class I Monoclonal Antibody, Unconjugated Clone W6/32 857.223.020

<b>Product Name</b>	Anti-Human HLA Class I Monoclonal Antibody, Unconjugated Clone W6/32
<b>Catalog Number</b>	857.223.020
<b>Description</b>	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC, IP
<b>Immunogen</b>	Membrane of human tonsil cells
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Stable at +2-8 °C for 12 months. DO NOT FREEZE.
<b>Note</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Isotype</b>	Mouse IgG2a Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
<b>Specificity</b>	Recognises an antigenic determinant shared among products of the HLA-A, -B and -C loci
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells