

Product Datasheet

Anti-Human CD123 Monoclonal Antibody, Azide Free Clone B-V8 857.230

Product Name	Anti-Human CD123 Monoclonal Antibody, Azide Free Clone B-V8
Catalog Number	857.230
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Recombinant human CD123
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile filtered through 0.22µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human CD123
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells