

Product Datasheet

Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-B53 857.260.005

Product Name	Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-B53
Catalog Number	857.260.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, WB
Immunogen	Recombinant human IL-1b
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human IL-1β
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells