

Product Datasheet

Anti-Human L1CAM Monoclonal Antibody, Unconjugated Clone B-L51 857.323.020

Product Name	Anti-Human L1CAM Monoclonal Antibody, Unconjugated Clone B-L51
Catalog Number	857.323.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human L1CAM/Fc
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human L1CAM
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells