

Product Datasheet

Anti-Porcine FSH Monoclonal Antibody, Azide Free Clone B-F52 857.620.005

Product Name	Anti-Porcine FSH Monoclonal Antibody, Azide Free Clone B-F52
Catalog Number	857.620.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Immunogen	Natural porcine FSH
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes both natural and recombinant porcine FSH
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells