

Product Datasheet

Anti-Human EpCAM Monoclonal Antibody, Azide Free Clone B-K46 857.625

Product Name	Anti-Human EpCAM Monoclonal Antibody, Azide Free Clone B-K46
Catalog Number	857.625
Description	Purified Monoclonal Antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant Human EpCAM/Fc
Alternative Names	CD326 TROP1
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises both natural and recombinant human EpCAM
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells