

Product Datasheet

Anti-Human CD146 Monoclonal Antibody, Azide Free Clone B-T46 857.630

| | |
|----------------------------|--|
| Product Name | Anti-Human CD146 Monoclonal Antibody, Azide Free Clone B-T46 |
| Catalog Number | 857.630 |
| Description | Purified Monoclonal Antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells |
| Species/Host | Mouse |
| Conjugation | Azide free |
| Tested Applications | WB, FC |
| Immunogen | Recombinant human CD146 |
| Alternative Names | MUC18 MCAM |
| Brand | Diaclone |
| Form/Appearance | Liquid |
| Concentration | 1.0 mg/ml (+/- 10 %) |
| Storage | Stable at +2-8 °C for 12 months. DO NOT FREEZE. |
| Isotype | Mouse IgG1 Kappa light chain |
| Clonality | Monoclonal |
| Buffer | Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. |
| Specificity | Recognizes both recombinant and natural human CD146 |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c node cells |