

Product Datasheet

Anti-Human IL-17C Monoclonal Antibody, Azide Free Clone B-Z17 857.710

Product Name	Anti-Human IL-17C Monoclonal Antibody, Azide Free Clone B-Z17
Catalog Number	857.710
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IL-17C
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises natural human IL-17C (but not all recombinants)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells