

Product Datasheet

Anti-Human CD63 Monoclonal Antibody, Azide Free Clone TS63 857.770

Product Name	Anti-Human CD63 Monoclonal Antibody, Azide Free Clone TS63
Catalog Number	857.770
Description	Purified Monoclonal Antibody. Source: Myeloma X63/AG8. x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	WB, FC, IP
Immunogen	Jurkat and HEL cell lines
Alternative Names	LAMP-3
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises the LIMP antigen, a 30-60 kDa (smear) protein
Hybridoma	Myeloma X63/AG8. x Balb/c spleen cells