

Product Datasheet

Anti-Human IL-5 Monoclonal Antibody, Azide Free Clone B-Z25 857.850.005

Product Name	Anti-Human IL-5 Monoclonal Antibody, Azide Free Clone B-Z25
Catalog Number	857.850.005
Description	Mouse monoclonal
Species/Host	Rat
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	COS-expressed recombinant hIL-5
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Rat IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Human IL-5