

Product Datasheet

Anti-Human TIM-3 Monoclonal Antibody, Unconjugated Clone B-K48 857.863.020

Product Name	Anti-Human TIM-3 Monoclonal Antibody, Unconjugated Clone B-K48
Catalog Number	857.863.020
Description	Purified Monoclonal Antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant Human TIM-3/Fc
Alternative Names	CD366 HAVCR2
Brand	Diaclone
Form/Appearance	Liquid
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognise natural (both soluble and membrane forms) and recombinant human TIM-3
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells