

Product Datasheet

Anti-Human ST2 Monoclonal Antibody, Azide Free Clone B-A52 857.870.005

Product Name	Anti-Human ST2 Monoclonal Antibody, Azide Free Clone B-A52
Catalog Number	857.870.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human ST2/Fc
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Store at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognize both natural soluble form and recombinant human ST2
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells