

Product Datasheet

Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-N42 857.910

Product Name	Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-N42
Catalog Number	857.910
Description	Purified Monoclonal Antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IFN alpha 2a
Alternative Names	IFNA2A Interferon alpha A
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes both natural and recombinant human interferon alpha 2a
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells