

Product Datasheet

Anti-HRP2 3902 SPTN-5 100560

Product Name	Anti-HRP2 3902 SPTN-5
Catalog Number	100560
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant HRP2 antigen, Medix Biochemica 610081
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.4
Cross Reactivity	N/D
Specificity	Antibody recognizes Plasmodium falciparum histidine-rich protein 2 (HRP2)
Shelf Life	36 months