

Product Datasheet

IgM Isotype Control , Unconjugated Clone PFR-03 857.913.020

Product Name	IgM Isotype Control , Unconjugated Clone PFR-03
Catalog Number	857.913.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Brand	Diaclone
Form/Appearance	Liquid
Storage	Store at +2-8 °C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgM
Clonality	Monoclonal
Buffer	Phosphate-buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Specificity	Reacts with undefined epitope on a plant pathogen