

Product Datasheet

Anti-Human HLA Class II (DR) Monoclonal Antibody, Azide Free Clone B-C10 857.930

Product Name	Anti-Human HLA Class II (DR) Monoclonal Antibody, Azide Free Clone B-C10
Catalog Number	857.930
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	RPMI 8866 and DAUDI cell lines
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the HLA Class II (DR) membrane
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells