

Product Datasheet

Anti-Human Galectin-9 Monoclonal Antibody, Azide Free Clone B-T33 857.940

Product Name	Anti-Human Galectin-9 Monoclonal Antibody, Azide Free Clone B-T33
Catalog Number	857.940
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC, ELISA
Immunogen	Recombinant human Galectin-9
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human Galectin-9 protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells