

Product Datasheet

Anti-h PCT 4003 SPTN-5 100562

Product Name	Anti-h PCT 4003 SPTN-5
Catalog Number	100562
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	CLIA, LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Anti-h PCT 4003 antibody recognizes the katalcalcin domain of procalcitonin, i.e. the region Asp(96) - Asn(116), as calculated from the peptide sequence of the mature PCT protein, without the signal sequence.
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.4-6.7
Cross Reactivity	Does not recognize canine procalcitonin.
Specificity	Antibody recognizes human procalcitonin (PCT)
Shelf Life	Unspecified