

## Product Datasheet

### Anti-h PCT 4003 SPTN-5 100562

<b>Product Name</b>	Anti-h PCT 4003 SPTN-5
<b>Catalog Number</b>	100562
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	LF, CLIA, FIA
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/- 10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	Anti-h PCT 4003 antibody recognizes the katalcalcin domain of procalcitonin, i.e. the region Asp(96) - Asn(116), as calculated from the peptide sequence of the mature PCT protein, without the signal sequence.
<b>Purity</b>	≥ 95 %
<b>Associated Products</b>	Recombinant PCT antigen, Medix Biochemica 610080
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.4-6.7
<b>Cross Reactivity</b>	Does not recognize canine procalcitonin.
<b>Specificity</b>	Antibody recognizes human procalcitonin (PCT)
<b>Shelf Life</b>	Unspecified