

Product Datasheet

Anti-Human ST2 Monoclonal Antibody, Azide Free Clone B-L42 857.950.005

Product Name	Anti-Human ST2 Monoclonal Antibody, Azide Free Clone B-L42
Catalog Number	857.950.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	WB, FC
Immunogen	Recombinant human ST2
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural membrane form and recombinant human ST2
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells