

Product Datasheet

Diaclone Lysing Buffer 858.040.125

Product Name	Diaclone Lysing Buffer
Catalog Number	858.040.125
Description	Erythrocyte lysing solution for preparation of leukocytes in whole blood for their analysis by flow cytometry.
Brand	Diaclone
Form/Appearance	Lyophilised
Storage	Store at room temperature. After reconstitution store at +2°C to 8°C.
Note	Reconstitute with 250 mL of distilled water, check and adjust the pH to 7.3. After reconstitution, the lysing solution is stable for 12 months.
Specificity	Erythrocyte lysing solution