

## Product Datasheet

### Diaclone Lysing Buffer 858.040.125

<b>Product Name</b>	Diaclone Lysing Buffer
<b>Catalog Number</b>	858.040.125
<b>Description</b>	Erythrocyte lysing solution for preparation of leukocytes in whole blood for their analysis by flow cytometry.
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Lyophilised
<b>Storage</b>	Store at room temperature. After reconstitution store at +2°C to 8°C.
<b>Note</b>	Reconstitute with 250 mL of distilled water, check and adjust the pH to 7.3. After reconstitution, the lysing solution is stable for 12 months.
<b>Specificity</b>	Erythrocyte lysing solution