

Product Datasheet

Murine IL-2 ELISA Kit 860.000.096

Product Name	Murine IL-2 ELISA Kit
Catalog Number	860.000.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	3.40 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine IL-2
Sample Types	Serum, Cell culture supernatant