

## Product Datasheet

### Murine TNF alpha ELISA Kit 860.040

<b>Product Name</b>	Murine TNF alpha ELISA Kit
<b>Catalog Number</b>	860.040
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	31.25-1000 pg/ml
<b>Sensitivity</b>	10.7 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	TNFa
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other murine cytokines
<b>Specificity</b>	Recognizes both natural and recombinant murine TNF-a
<b>Sample Types</b>	Serum, Cell culture supernatant