

Product Datasheet

Murine IFN gamma ELISA Kit 860.050

| | |
|---------------------|--|
| Product Name | Murine IFN gamma ELISA Kit |
| Catalog Number | 860.050 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 31.25-1000 pg/ml |
| Sensitivity | 15 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | IFN-g, IFN-gamma |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other murine cytokines |
| Specificity | Recognizes both natural and recombinant murine IFN-g |
| Sample Types | Serum, Cell culture supernatant |