

# Product Datasheet

## Murine IFN gamma ELISA Kit 860.050.048

Product Name	Murine IFN gamma ELISA Kit
Catalog Number	860.050.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine IFN-g
Sample Types	Serum, Cell culture supernatant