

# Product Datasheet

## Murine GM-CSF ELISA Kit 860.060.048

Product Name	Murine GM-CSF ELISA Kit
Catalog Number	860.060.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	9 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine GM-CSF
Sample Types	Serum, Cell culture supernatant