

Product Datasheet

Murine CD138 ELISA Kit 860.090

Product Name	Murine CD138 ELISA Kit
Catalog Number	860.090
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.5-16 ng/ml
Sensitivity	0.35 ng/ml
Tested Applications	ELISA
Alternative Names	Syndecan-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other molecules
Specificity	Recognizes recombinant and natural murine soluble Syndecan-1
Sample Types	Cell culture supernatant, serum, plasma, buffered solution